



STATEWIDE

HOME HEALTH CARE




WHEELCHAIRS

Your home. Your health. We care.

LIGHTWEIGHT WHEELCHAIRS

The S-ERGO series of Lightweight wheelchairs feature Karma's unique S-Ergo shaped seating system.

The ergonomic 'S' shaped seat helps spread the load and prevents the user slipping down in the chair. The seat and backrest are made from a mesh upholstery that is great for keeping cool on summer days. The S-Ergo seat has been pressure-mapped for optimum shape. In addition the seating system benefits from the unique anti-bacterial technology built into the upholstery that provides protection against odour, staining and deterioration caused by bacteria, fungus and other microorganisms.



150834B

150835

Scan here to learn more!



ERGO LITE RANGE

The Ergo Lite Deluxe is the lightest wheelchair in our range. Equipped with swing away removable leg rests and double cross bracing for strength and durability.

KARMA ERGO LITE MODELS		
Self Propelled Models	41 cm x 41 cm (16" x 16")	46cm x 41 cm (18" x 16")
Seat Height	48 cm	48 cm
Product Weight	10.2 kg	10.3 kg
Max User Weight	100 kg	100 kg
Transit Models	41cm x 41cm (16" x 16")	46cm x 41cm (18" x 16")
Seat Height	48cm	48 cm
Product Weight	8.7 kg	8.8 kg
Max User Weight	100 kg	100 kg



150813

150814 TRANSIT

Scan here to learn more!



S-ERGO 125

In addition to the S-Ergo seating system, the S-Ergo 125 also has seat and backrest upholstery and pads treated with Aegis Microbe shield, to reduce stains and odour.

The streamline oval tubes not only make the wheelchair look great, but also add to the strength of the wheelchair structure. Complete with swing away detachable footplates.

KARMA S ERGO 125 LITE MODELS			
Self Propelled Models	41cm x 43cm (16" x 17")	46cm x 43cm (18" x 17")	51cm x 43cm (20" x 17")
Seat Height	51 cm	51cm	51 cm
Product Weight	14.2 kg	14.5kg	14.8 kg
Max User Weight	115 kg	115kg	115 kg
Transit Models	41cm x 43cm (16" x 17")	46cm X 43cm (18" x 17")	51cm x 43cm (20" x 17")
Seat Height	50 cm	50 cm	50 cm
Product Weight	13.1 kg	13.1 kg	13.1 kg
Max User Weight	115 kg	115 kg	115 kg

WHEELCHAIRS

Our range of wheelchairs cater to a wide range of needs, offering versatility for both the user and carers.



15-0869SP



15-0870T



15-0872M



15-0872OG

VIDA FOLDING WHEELCHAIR

Super lightweight transport wheelchair ideal for quick trips. Ergonomic backrest and seat shape follow the contours of the user to maximise comfort and postural support.

Attendant operated hand brakes with parking lock. Removable leg rests and folding back for compact storage. Rear storage pocket for essentials.

Colours available:

- Mustard – 40 cm seat width
- Magenta – 45 cm seat width
- Teal – 45 cm seat width
- Olive – 50 cm seat width

Available in Self Propelled and Transit Models.

SELF-PROPELLED	
Seat Width	40 / 45 / 50 cm
Product Weight	11.7 kg
Max User Weight	135 kg

TRANSIT	
Seat Width	40 / 45 / 50 cm
Product Weight	9.6 kg
Max User Weight	135 kg



15-0846a



15-0804

SOCIALITE WHEELCHAIR

Ultra light weight transport chair perfect for short trips. Padded and breathable seat and back provide optimum comfort and postural support. Attendant operated hand brakes with parking lock. Folds for easy transportation.

Available in Self Propelled model.
Colours available Red or Silver.

Seat Width	45 cm
Product Weight	8.1 kg
Max User Weight	115 kg

EVOKE WHEELCHAIR

Designed for the user who wants a lightweight, highly adjustable manual wheelchair. The contoured padded overlay offers the user increased comfort and is ideal for those who are in their wheelchair for longer periods of time. Swing-in/out, height adjustable and removable foot rests allows optimal positioning for the changing needs of the user. Crash Tested.

Available in 45 cm and 50 cm seat widths.

Seat Width	45 cm	50 cm
Product Weight	15 kg	15 kg
Max User Weight	140 kg	140 kg



15-0806SP

DAYS WHIRL

Designed to be a light, practical and easy to transport wheelchair, making it an ideal occasional use chair for both indoor and outdoor use. Light construction and ease to fold for storage or transport make it a very versatile chair.

Seat Height	50 cm
Seat Width	45 cm
Product Weight	16.8 kg
Max User Weight	120 kg

WHEELCHAIRS

HORIZON ALUMINIUM WHEELCHAIR

Lightweight wheelchair suitable for day trips.

Standard Features

Quick release wheels

Lift up armrests

Attendant handbrakes

User lap belt and anti tippers



15-0806A



15-0837

Available in Self Propelled and Transit Models.

SELF-PROPELLED	
Seat Width	46 cm
Product Weight	15 kg
Max User Weight	125 kg

TRANSIT	
Seat Width	46 cm
Product Weight	12 kg
Max User Weight	125 kg

OTTOBOCK START M2 SERIES

The Otto Bock Start M2S is the entry-level model of the Start series. The chair is equipped with split, pivoted-away footrests. In addition, the side panels are equipped with height-adjustable armrests, so you are always comfortably supported.

The Otto Bock Start can be adapted effortlessly to individual body size and needs, allowing the user to be optimally positioned and supported at all times. The back height is 42.5 inches set fixed. 20 variants of the model are ready for delivery, the selection is made with a simplified order form. You just have to indicate the appropriate seat width and choose a variant. The integrated tip assists the tilting of the wheelchair, making it easier to negotiate stairs.

Wheelchair Features Include:

- Increased weight capacity to 130kg + Crash Tested
- Horizontal Rear Wheel Adjustment on all models
- Patented Quick Depth Adjustment (+/- 1 inch)
- Anti-Bacterial Push Handles (ProvaMed®)
- Improved Cross Brace Clips
- Asymmetric, reversible footplates

ottobock.
A new generation
of wheelchairs



15-0854

**Weight from 14.7 kg
(Transport weight from 9.6 kg)**

Available in Three Models

Standard, Adjustable and Configurable.

Seat Widths Available

Standard	405/455/50 mm
Adjustable	380/405/430/455/480 mm
Configurable	380/405/430/455/480/505 mm

Options

Side Panel with Height Arm Adjustment
One Arm drive
Armrest with swivel unit
Elevating Footrest
Amputating leg rest
Add on Drives (i.e the new support)

POWER ASSIST

The Power Assist attaches directly on to most manual wheelchairs. It features an easy-to-use squeeze trigger mechanism, along with a speed dial, enabling users to customise speed settings.

The Power Assist is designed to work in conjunction with most wheelchairs with seat widths of 16" – 20". (wheelchair and user combined).

- The Power Assist has a maximum speed of up to 5km/h (3.1 mph).

Maximum User Weight 120 kg

**Wheelchair not included.*



19-0122

STANDARD WHEELCHAIRS

BREEZY® BASIX 2

This wheelchair represents improved adjustability and rigidity, combined with a modern, clean look. Improvements include greater seat depth adjustment (41 cm to 46 cm), a tension adjustable backrest which is extendable from 41 cm to 46 cm, and a more robust, rigid castor attachment (0, 3 and 6 degrees). Highly efficient and easy to use. We also hire these wheelchairs in 14 and 16 for Children and 18 inch for Adults.

Standard Features

- Transit Model Available
- Castor Fork Adjustment
- Axle Plate Adjustment
- Height Adjustable handles
- Elevating Leg Rest
- Rounded Arm Rests

15-0857



Available in the following seat widths

Seat Width	Weight Capacity
35.6 cm (14 Inch)	125 kg
40.64 cm (16 Inch)	125 kg
45.72 cm (18 Inch)	125 kg
50.8 cm (20 Inch)	125 kg

BREEZY® EC SERIES WHEELCHAIRS

The dependable, industry-standard Breezy® EC wheelchairs have your mobility needs covered.

Standard Features

- Swing-away foot plates
- Dual axle allows the seat to be lowered to hemi height
- Chrome-plated, carbon steel frame is built to last
- Wheels have moulded contours to add strength and durability
- Full size chart pocket
- Industry-compatible front riggings

Available in the following seat widths

Seat Width	Weight Capacity
45.72 cm (18 Inch)	113.5 kg
50.8 cm (20 Inch)	135 kg
55.84 cm (22 Inch)	205 kg
60.96 cm (24 Inch)	205 kg



15-0802

TRANSIT BUDDY

The Transit Buddy is an excellent lightweight travel chair, designed for convenience and ease of use. Weighing just 8.4kg

Standard Features

- Foldable footrest
- No rust aluminium frame
- No assembly required
- Lockable handbrake system

Seat Width	15"
Product Weight	4.8 kg
Max. User Weight	100 kg




15-0802

ACCESSORIES

Leg Extensions

Left - 150808L
Right - 150808R



Stump Support

Left - 15080L
Right - 15080R



Seat Belt Assembly

150913



Oxygen Bottler Holder

Suitable for Breezy EC Wheelchair
150901



TILT-IN-SPACE WHEELCHAIRS



15-0842

REHAB RS CLASSIC TILT-IN-SPACE WHEELCHAIR

The Rehab RS is designed for a user who has limited ability to reposition themselves in a device, and as a result is at risk of pressure injuries.

- Tilt (26°) and backrest recline (30°) facilitates weight re-distribution and the preservation of skin integrity. It also optimises user comfort, decreases user fatigue and assists with head control
- Adjustable seat depth and width
- Height adjustable and removable armrests
- Adjustable headrest
- Heavy duty, padded, elevating leg rests (std)
- Crash Tested – to AS/NZ3616.196

Maximum User Weight 180 kg

BREEZY RELAX TILT-IN-SPACE WHEELCHAIR

The Breezy Relax is a multi-function tilt and recline (including anterior tilt) wheelchair with a wide range of features and adjustments.

- Simply adjusted without special tools
- Comes standard with anatomically shaped seat and backrest
- Adjustable seat height and depth, adjustable height removable armrests with sliding arm pads
- Elevating detachable leg rests with angle adjustable footplates, reclining height adjustable backrest with removable adjustable three part headrest
- Available in transit or self propel models

Maximum User Weight 140 kg

Scan here to learn more!



15-0849

BREEZY CIRRUS TILT IN SPACE WHEELCHAIR

The Cirrus G5 provides high levels of comfort combined with an infinite level of adjustment possibilities including comfortable upholstery and arm rests that follow the back when reclining offering optimal support and comfort all day long.

- Numerous adjustment options that ensure a good and stable seat position
- The back recline and tilt provides an easy and comfortable way to change position throughout the day
- Adjustments can be made while the user is sitting in the wheelchair, and they do not require the use of any tools.

Maximum User Weight 135 kg



15-0861



Scan here to learn more!

ROLLZ MOTION ELECTRIC

Whether you need accompaniment or prefer to go solo, you can rely on the Rollz Motion Electric. This energy-efficient option is for those looking to maintain their mobility.



The Rollz Motion Electric seamlessly switches between four mobility-aid options:

1/ Rollator

Walk Safely and with confidence.

Your safety is a top priority, which is why the Rollz Motion Electric has been designed with a sturdy and reliable frame with a stylish and modern design that looks great wherever you go.

2/ Transport chair

The easy-to-use conversion seat and footrests means you can transform your rollator into a wheelchair in seconds, without any additional tools or assistance.

3/ Power Assist

When configured as a transport chair, The Rollz Motion Electric features a power assist that provides an electric push support to ease pushing uphill on ramps, slopes, and hills. This is especially beneficial for the person assisting the user, as it reduces the physical effort required to push the wheelchair.

4/ Powerchair

The Rollz Motion Electric is also equipped with a joystick controller that allows for independent use of the wheelchair. By pressing the On/Off button, the joystick is activated and users can control the device's movement and speed. The +/- buttons allow for a change in gears and the horn button can be used to alert others of your approach or presence.




Scan here
to learn
more!

Rollator Weight	19.6 kg
Battery Weight	3 kg
Maximum Load	100 kg
User Height	1.6 - 1.9 m
Seat Height	55 cm
Handle Height Range	88 - 99 cm
Overall Width X Length	74 x 71 cm

- Push support option to ease the assisted transport chair
- Joystick controller for independent driving in the electric wheelchair
- 15 km range
- Maximum speed of 6 km/h
- Easy foldability that allows compact transportation
- Low-vibration air tyres

Learn more about other popular products that may support your home health care journey.



Statewide Home Health Care
909 Nepean Highway, Bentleigh 3204
T. (03) 9591 6234 | sales@shhc.com.au | shhc.com.au



Hire and sales also available on the Mornington Peninsula:

Peninsula Home Health Care
21 Sir Laurence Drive, Seaford 3198
T. 9786 7004 | sales@phhc.com.au | phhc.com.au

 @shhc_cares  @Statewide Home Health Care  @StatewideHomeHealthCare