



# STATEWIDE

HOME HEALTH CARE



## POWERED MOBILITY

Your home. Your health. We care.

# TRANSPORTABLE MOBILITY SCOOTERS

## PRIDE GO-GO BOOT SCOOTER

The Boot Scooter has one-hand feather-touch disassembly and a convenient drop-in battery box that makes transport and travel worry-free. Disassemble into 5 pieces for easy transportation and storage.

|                                |           |
|--------------------------------|-----------|
| Max Weight Capacity            | 118 kg    |
| Max Speed                      | 6.4 km/h  |
| Travel Distance On Full Charge | 11.5 km** |



19-0401



19-0114 LX4

Scan here to learn more!



## GO-GO LX WITH CTS

The Go-Go® LX is the most advanced travel scooter available.

- Independent front and rear Comfort-Trac Suspension (CTS)
- CTS offers a new level of performance giving a smooth and responsive ride
- Disassembles into 5 pieces for easy transportation and storage

|                                |           |
|--------------------------------|-----------|
| Max Weight Capacity            | 136 kg    |
| Max Speed                      | 6.4 km/h  |
| Travel Distance On Full Charge | 16.9 km** |

## YOGA FOLDING

Fold and Roll, or store in a suitcase. Ideal for both indoor and outdoor use.

- High ground clearance great for cruises
- Height retractable and angle adjustable tiller
- Strong light-weight aluminium frame designed for easy folding and transporting
- Supplied with a strong sturdy hard case to keep the scooter safe when travelling

Scan here to learn more!



\*\* Depends on user weight, terrain type, battery charge, battery condition and tyre condition



# TRANSPORTABLE MOBILITY SCOOTERS



190113

## GO-GO ULTRA X

The Ultra X has a compact, lightweight design that ensures you always stay independent.

- Disassembles into 5 easy pieces for convenience
- Available in 3 and 4 wheel model
- Black solid non scuffing tyres
- Lightweight 12amp battery box
- Complete with 12amp battery

|                                |          |
|--------------------------------|----------|
| Max Weight Capacity            | 118kg    |
| Max Speed                      | 6km/h    |
| Travel Distance On Full Charge | 12.8km** |

**PRIDE**  
MOBILITY  
LIVE YOUR BEST\*



## SONIC SCOOTER

This portable and durable scooter quickly dismantles into five lightweight pieces, no tools required.

- The Sonic's compact size is the perfect choice for your on-the-go lifestyle
- Easy transport and storage
- Puncture-proof tyres and a detachable battery pack
- Adjustable angle tiller with delta handle

|                                |          |
|--------------------------------|----------|
| Max Weight Capacity            | 136 kg   |
| Max Speed                      | 6.8 km/h |
| Travel Distance On Full Charge | 20 km**  |



19-0400

## THE ATTO

A versatile mobility scooter ideal for use on public transport, including ship and air travel. The "ATTO" can be easily parked and folded for transport.

- Assembled Weight 30kg. When Split, each half weighs 15 kg

|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 100 kg  |
| Max Speed                      | 6 km/h  |
| Travel Distance On Full Charge | 16 km** |

Scan here to learn more!



Scan here to learn more!



1190300



\*\* Depends on user weight, terrain type, battery charge, battery condition and tyre condition

## MEDIUM SCOOTERS

### ECO 745



Stylish, eye catching design with attractive modern colour range.

- Available in green, yellow or red
- Loop handlebars offer comfortable finger or thumb throttle operation
- Adjustable high back contoured seat with simple operation of slide and swivel mechanisms to allow maximum user comfort
- Adjustable power saving LED headlight, indicators front and rear, rear tail lights and reverse warning
- 2 Pole Scooters – suitable for flat surfaces and lower speeds  
Ideal for indoor and occasional use

|                                |        |
|--------------------------------|--------|
| Max Weight Capacity            | 136kg  |
| Max Speed                      | 10km/h |
| Travel Distance On Full Charge | 20km** |



Scan here to learn more!



19-0220

### FENDE S946

The Merits Fende will negotiate the footpaths with ease and is still compact enough to take you into the shopping centre

- Leather feel contoured high back seat with slide, swivel, recline and adjustable armrests
- Unique battery storage allows for greater floor space and a safe secure lockable storage area for personal items
- New modern LCD dashboard display offers clearer, more easy to read battery gauge, speed, time & temperature
- Built in speed reduction switch which is activated, when turning offering exceptional safety
- 2 Pole Scooters – suitable for flat surfaces and lower speeds  
Ideal for indoor and occasional use



19-0221

|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 150 kg  |
| Max Speed                      | 10 km/h |
| Travel Distance On Full Charge | 30 km** |



Scan here to learn more!

### MERITS EXPLORER

Streamline design with a sporty feel, available in red, silver or kite blue.

- Full suspension system with shock absorbers that delivers a smooth, comfortable ride
- Moder LED lighting system, Digital LCD dashboard: time display, battery gauge, temperature, maintenance messages, speed display and odometer
- Extra storage space beneath tiller
- Half speed safety switch which operates automatically when turning to provide extra safety
- 4 Pole Scooter – suitable for inclines, outdoor and more frequent use

|                                |            |
|--------------------------------|------------|
| Max Weight Capacity            | 159 kg     |
| Max Speed                      | 12 km/h    |
| Travel Distance On Full Charge | 40-50 km** |



19-0223

\*\* Depends on user weight, terrain type, battery charge, battery condition and tyre condition

## MEDIUM SCOOTERS



### 745 PLUS

This new design offers generous foot space & leg room and large battery capacity.

- Available in blue, silver or yellow
- Adjustable high back contoured seat with simple operation of slide & swivel mechanism to allow maximum user comfort
- Adjustable power saving LED headlight, front and rear indicators, rear tail and brake lights
- Built in speed reduction switch, activated when turning, offering exceptional safety
- 4 Pole Scooter – suitable for inclines, outdoor and more frequent use



19-0222

|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 150 kg  |
| Max Speed                      | 10 km/h |
| Travel Distance On Full Charge | 30 km** |

### PRIDE 130XL SCOOTER



Equipped with a robust 4-pole motor, this scooter is suitable for both paved paths and rugged surfaces with front and rear suspension provide a stable, comfortable ride over outdoor terrain.

- Large, tiller-mounted front basket, high-back reclining seat and rear view mirror
- Large pneumatic tyres for excellent outdoor performance
- Bright, long-lasting, front and rear LED lighting and directional indicators
- Adjustable, delta tiller with ergonomic wraparound handles letting you operate the scooter with one, easily accessible tie down points
- 4 poles
- Available in Red and Silver



|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 158 kg  |
| Max Speed                      | 10 km/h |
| Travel Distance On Full Charge | 25 km** |

# LARGE SCOOTERS

## INTERCEPTOR S2

- Streamline design with a sporty feel
- Fully adjustable suspension system that delivers a smooth and comfortable ride
- Modern LED lighting system delivers brilliant light with less power consumption
- Digital LCD dashboard display including a time display, battery gauge, temperature, maintenance messages, speed display and odometer
- Extra storage space beneath tiller
- 4 Pole Scooter – suitable for inclines, outdoor and more frequent use

|                                |              |
|--------------------------------|--------------|
| Max Weight Capacity            | 200 kg       |
| Max Speed                      | 8.9 km/h     |
| Travel Distance On Full Charge | 40 - 60 km** |

19-0304



Scan here to learn more!



## PRIDE 140XL

Equipped with a robust 24- volt DC motor, full suspension, hydraulic-sealed brake system, and 13” pneumatic tyres, this scooter is suitable for any terrain.

- Large, tiller-mounted front basket and High Back Pillo top Seat
- Front and rear suspension and sealed hydraulic brakes
- Large 13” pneumatic tyres for excellent outdoor performance
- Bright, long-lasting, front and rear LED lights and directional indicators
- Easily accessible tie down points (for the transport of unoccupied scooter) front and rear LED lights and directional indicators and rearview mirror
- Turning Radius 210cm, 13cm Ground Clearance
- Available in Red, Yellow or Silver

|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 181 kg  |
| Max Speed                      | 10km/h  |
| Travel Distance On Full Charge | 35 km** |

190108

**PRIDE**  
MOBILITY  
LIVE YOUR BEST®



## ACCESSORIES

|                                                                                    |                                                                                     |                                                                                     |                                                                                      |                                                                                       |
|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |  |  |
| <b>Rear Pod</b><br>190211                                                          | <b>Cup Holder</b><br>190214                                                         | <b>Walking Frame Holder</b><br>190204                                               | <b>Rear Bag</b><br>190206                                                            | <b>Safety Flag</b><br>190202                                                          |



## POWER CHAIRS

### KARMA KP-25.2

The KP-25.2 is a foldable electric wheelchair making transport convenient.

- Built-in multi-functional controller
- Flip-back and height adjustable armrests
- Fire retardant upholstery

|                                |            |
|--------------------------------|------------|
| Seat Dimensions                | 46 x 46 cm |
| Seat Height                    | 50 cm      |
| Max Weight Capacity            | 135 kg     |
| Max Speed                      | 10 km/h    |
| Travel Distance On Full Charge | 30 km**    |



Scan here to  
learn more!

19-0225



### GO CHAIR

The first travel chair that is super portable, with manoeuvrability and style at a low cost.

- Conveniently disassembles into 5 pieces
- Compact design for tight turning radius
- Indoor and smooth surface use only

|                                |            |
|--------------------------------|------------|
| Max Weight Capacity            | 113 kg     |
| Max Speed                      | 5.65 km/h  |
| Travel Distance On Full Charge | 12.87 km** |



19-0103



### FOLD & GO POWER WHEELCHAIR

Folds in one easy step, lightweight.

- Allows manoeuvrability in tight spaces.
- Adjustable footplate and backrest.

|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 114kg   |
| Max Speed                      | 4.7km/h |
| Travel Distance On Full Charge | 11km**  |



Scan here to  
learn more!

19-0226



\*\* Actual distance in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

## POWER CHAIRS - ULTRA LIGHT



19-0402

### MERITS LITE RYDER

- Super Lightweight carbon fibre frame, folds in one easy step for transportation
- Weight w/out battery 11.85 kg
- Complete Overall Weight with battery 13 kg
- Comfortable padded seat & backrest
- USB fitting in Joystick & in removable battery pack

|                                |          |
|--------------------------------|----------|
| Max Weight Capacity            | 100 kg   |
| Max Speed                      | 4.5 km/h |
| Travel Distance On Full Charge | 10 km**  |

**merits**  
**AUSTRALIA**



19-0405

### QUICKIE Q50 R CARBON FIBRE POWER CHAIR

The Q50 R Carbon is the perfect travelling partner. Weighing a mere 14.5 kg (without batteries) and fabricated from premium carbon fibre, making this an effortless powerchair to lift.

- Foldable in seconds
- Compact – 37 cm x 57 cm x 76 cm
- Approved for Air Travel

|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 136 kg  |
| Max Speed                      | 6 km/h  |
| Travel Distance On Full Charge | 24 km** |

Scan here to  
learn more!



**SUNRISE**  
**MEDICAL**



19-0404

### FLYTEMAX2 POWERCHAIR

- Sleek High Strength Carbon Fibre Frame
- Fixed or Remote Attendant Controller
- Weight without Batteries 15.6k
- Complete Overall weight with batteries 17.4 kg
- Easy to manoeuvre
- Ergonomic Seat cushion for superior comfort
- Intuitive controller knob for seamless navigation
- Sturdy footrest for optimal support
- Seat width 500 cm

|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 150 kg  |
| Max Speed                      | 6 km/h  |
| Travel Distance On Full Charge | 20 km** |

**TopGun**  
**Mobility**  
Mobility & Homecare Specialists



\*\* Actual distance in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.



# MAVERICK POWER CHAIR RANGE

The Merits Maverick range mid-drive, 6 wheel design delivers both manoeuvrability in tight indoor spaces and stability outdoors.



## All Merits Maverick models feature

- Dual in-line motors offering incredible efficiency
- Low centre of gravity provides a smooth and stable ride even over rough terrains
- Swing away joystick mount
- Height-adjustable seat with height and width adjustable armrests
- A number of additional seating options available
- Optional left-handed or Attendant Control

## MAVERICK 10 P326A

- Indoor and outdoor use for the small to medium user
- Range of 28km\*

### 136 kg user capacity

|                                |        |
|--------------------------------|--------|
| Max Weight Capacity            | 136kg  |
| Max Speed                      | 8km/h  |
| Travel Distance On Full Charge | 28km** |



19-0302-10

## MAVERICK 12 P327

- Indoor and outdoor use for the medium to larger user
- Range of 32 km\*

### 205 kg user capacity

|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 205 kg  |
| Max Speed                      | 8 km/h  |
| Travel Distance On Full Charge | 32 km** |
| Turn Radius                    | 53 cm   |
| Ground Clearance               | 8 cm    |



15-0302-12

## MAVERICK 14 P330

- Indoor and outdoor use for the medium to larger user
- 4 Pole Technology motors for additional outdoor use
- Range of 40 km\*

### 205 kg user capacity


|                                |         |
|--------------------------------|---------|
| Max Weight Capacity            | 205 kg  |
| Max Speed                      | 7 km/h  |
| Travel Distance On Full Charge | 40 km** |
| Turn Radius                    | 53 cm   |
| Ground Clearance               | 5.6 cm  |



15-0302-14

\*\* Actual distance in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement.

Learn more about other popular products that may support your home health care journey.



**Statewide Home Health Care**  
909 Nepean Highway, Bentleigh 3204  
T. (03) 9591 6234 | [sales@shhc.com.au](mailto:sales@shhc.com.au) | [shhc.com.au](http://shhc.com.au)



Hire and sales also available on the Mornington Peninsula:

**Peninsula Home Health Care**  
21 Sir Laurence Drive, Seaford 3198  
T. 9786 7004 | [sales@phhc.com.au](mailto:sales@phhc.com.au) | [phhc.com.au](http://phhc.com.au)

 @shhc\_cares  @Statewide Home Health Care  @StatewideHomeHealthCare