



PATIENT TRANSFER

Your home. Your health. We care.

LIFTS, SLINGS & MANUAL AIDS

A range of easy to operate mobile lifters that allow patients to be lifted and transported from the floor, bed or chair with ease and safety.

LIFTING HOIST – ALTO LIFT MKII 200

The ALTO LIFT MKII 200 is a medium sized heavy duty mobile patient lifting hoist with a safe working load of 200 kg and a compact ergonomic design.

This lifter is ideally suited for institutional use or as a compact heavy duty home care hoist for the bariatric patient.

This ergonomic lifting is constructed with a powder coated aluminium frame for ease of manoeuvrability and easy to wipe down for cleaning.

Safe Working Load	200 kg
Base Length	1150 cm
External Base Width Closed	700 mm
External Base Width Open	1130 mm
Internal Base Width Closed	600 mm
Internal Base Width Open	1030 mm
Base Width Adjustment	Electric
Under Bed Clearance	130 mm
Minimum Height Yoke	300 mm
Maximum Height Yoke	1920 mm
Average Lift Range*	1620 mm
Lifter Weight	46 kg
Minimum Storage Height	1350 mm
Under Base Height	55 mm
Maximum Hoist Reach	1060 mm
Minimum Hoist Reach	610 mm
Hoist Reach At Maximum Height	515 mm



21-0502

INVACARE BIRDIE COMPACT 150

Designed to fold easily and without the need for tools, the Birdie Compact can be stored or transported comfortably when space is limited.

The compact base and large castors manoeuvre smoothly around beds, in bathroom areas and through narrow doorways. The Birdie Compact is ideal for everyday lifting in the home or an aged care facility environment. Compact without compromising on comfort.

Safe Working Load	150 kg
Total Width Closed (External Measure)	55 cm
Total Length	100 cm
Total Internal Width (Base Closed/Open)	46 cm/89 cm
Lifting Range (From Sling Attach Point)	53 cm-160 cm
Turning Radius	170 cm
Base Height (Top Of Legs)	11.5 cm
Total Weight (Inc. Spreader Bar)	36 kg
Weight Of Heaviest Part (Base)	16.5 kg
Manual Emergency Lowering	No
Electric Emergency Lowering	Yes (up and down)
Castors	Dual 100 mm
Charging	Internal and External Std
Overall Height	112 cm

**Folds away
neatly for
easy storage**



210101B

LIFTS, SLINGS & MANUAL AIDS

FORTE 320 HEAVY DUTY LIFTING HOIST

With a safe working load of 320 kg, this extremely heavy duty mobile patient lifting hoist provides full bariatric capability whilst remaining a compact and ergonomic design. Manoeuvrability and exceptional under bed clearance are key design features.

FORTE 320 provides superior patient lift space through the combination of the recessed actuator and the use of SPACE LINK™ technology. The aesthetic design of the concealed electric leg mechanism provides infection control advantages.

FORTE 320 is the lifter of choice for institutional use where a high capacity bariatric lifter is required.

Safe Working Load	320 kg
Base Length	1433 mm
External Base Width Closed	715 mm
External Base Width Open	1205 mm
Internal Base Width Closed	590 mm
Internal Base Width Open	1080 mm
Base Width Adjustment	Electric
Under Bed Clearance	126 mm
Minimum Yoke Height	260 mm
Maximum Yoke Height	1955 mm

21-0503



INVACARE ISA XPLUS STANDING HOIST

The ISA XPlus is able to take on a user weight of 200 kg and is designed to suit a variety of users who require assistance during transfers, but can still weight bear. It provides safety, comfort and a more pleasant experience for clients.

The ISA XPlus features ergonomic hand grips that fit the natural shape of the hand and position of the hand and position of the wrist. Both the lifting arm and lower leg support can be adjusted to different height settings to suit the user's needs.

Other features include a push bar to ensure the smoothness and ease of mobility as well as a non-slip and large foot pedal for enhanced caregiver comfort.

- Ergonomic hand grips - its unique design fits perfectly with the natural shape of the hand for a safer grip
- Extendable lifting arm - The TelescoLift can be adjusted to 9 different settings to suit individual's height, thus offering maximum flexibility
- Lower leg support - Providing 6 height positions along with a swivel ErgoSupport option for optimum comfort
- Foot plate - 3 heights available with removable heel block and round shape for toilet approach
- Push bar and foot pedal
- Sling hooks
- Leg band
- Curved legs

Length	1280 mm
Min Base Width	640 mm
Inside Base Width	1020 mm
Lifting Range	505-675 mm
Castor Size	100,125 mm
Ground Clearance	35 mm
Max User Weight	200 kg
Unit Weight	52 kg
Battery Time	82 cycles
Overall Height	112 cm



210201B

LIFTS, SLINGS & MANUAL AIDS

SALSA STANDING HOIST

The SALSA mobile standing hoist is designed to assist individuals with limited mobility in transferring from a seated to a supported standing position. This electric-operated hoist features ergonomic design elements that prioritise both patient and carer comfort. With its small footprint and low weight, the Salsa 200 is highly manoeuvrable.

Safe Working Load	200 kg	Minimum Height Lift Arm	880 mm
Base Length	1010 mm	Maximum Height Lift Arm	1685 mm
External Base Width Closed	665 mm	Under Base Height	55 mm
External Base Width Open	1030 mm	Lifter Weight	51 kg
Internal Base Width Closed	565 mm	Minimum Storage Height	1190 mm
Internal Base Width Open	930 mm	Footplate	Removable
Base Width Adjustment	Electric	Kneepad Height	Adjustable
Under Bed Clearance	130 mm		



Scan here to learn more!



21-0106

ETAC TURNER

A practical aid which supports and activates the patient when transferring from one seated position to another. It offers security for the person being transferred and will not strain the caregivers' back.

Easy and safe transfers

- The Turner reduces the risk of falling, and it also allows the caregiver to work in an ergonomically correct position
- Unique patient turner and transfer aid that facilitates seat to seat transfers, for patients who are able to bear their own weight when standing
- Allows carers to work ergonomically, without strain to back or arms while still maintaining eye contact
- Promotes continued mobility as patient actively moves to a standing position
- Movement up and down controlled by the patient for enhanced security
- wLow footplate (13 mm high) makes foot placement easy while non-slip surface offers additional security



21-0402

OXFORD TRANSFER AID

The Oxford Up 200 is a uniquely designed manual stand-aid and transfer device, suitable for those who require assistance when standing but are able to contribute effort to the process.

By encouraging user participation, the Oxford Up 200 can play a supportive role in patient rehabilitation, function, and independence.

Safe Working Load	200 kg (31 st)
Overall Height (A)	1093 mm
Overall Length (B)	960 mm
Overall Width (C)	633 mm
Legs Open - External Width (D)	840 mm
Legs Open - Internal Width	720 mm

MANUAL HANDLING

Our range of manual handling accessories can be an invaluable tool for any domestic care situation, designed with the carer in mind.

A/ Hygiene sling with head support

The Hygiene with Head Support is specifically designed for toilet transfers for patients with less control, requiring additional head support. For Yoke frame lifters only. The open area from mid back to mid-thigh for hygiene access, buckle waist strap for added security and increased infection control, generous padding under arms.

B/ General Purpose Sling

Available in four sizes in a choice of fabric or mesh.

The General Purpose sling provides easy fitting and good support for patients with some upper body control.

C/ Paediatric Sling

Available in three sizes in mesh only. The sling provides easy fitting and good support for children or small adults. Designed for general lifting with children or small adults.

D/ Standup Slings

Standup Slings are designed to lift the patient comfortably from a seated position into a relaxed upright position.

The high-cut, padded Stand Up Sling supports the patient from waist to shoulder while allowing the arms to move freely so that the hands can reach naturally to grip the support rails as the patient rises.

E/ General Purpose Sling with Head Support

Available in five sizes in a choice of fabric or mesh. The sling provides easy removal of clothing or toilet transfers.


Recommended for general toilet transfers for the less secure patient.



General Purpose

	XS	S	M	L	XL	XXL
A	790	790	890	1110	1240	1400
B	210	230	240	260	270	290
C	470	460	530	650	730	750
D	1120	1130	1190	1390	1520	1590
FABRIC	GP01XS	GP01S	GP01M	GP01L	GP01XL	GP01XXL
MESH	GP02XS	GP02S	GP02M	GP02L	GP02XL	GP02XXL


YOKE ATTACHMENT (SPREADER BAR)



With head support

	XS	S	M	L	XL	XXL
A	820	820	920	1110	1270	1360
B	210	230	240	260	270	290
C	750	780	920	1000	1020	1050
D	1310	1380	1410	1630	1680	1780
FABRIC	GPHS01XS	GPHS01S	GPHS01M	GPHS01L	GPHS01XL	GPHS01XXL
MESH	GPHS02XS	GPHS02S	GPHS02M	GPHS02L	GPHS02XL	GPHS02XXL

YOKE ATTACHMENT (SPREADER BAR)



MANUAL HANDLING

ETAC SLING ON

Primarily designed for easy insertion of a hoist sling in bed without the need for rolling or moving the patient.

The Etac Sling On is an optimal tool for users that need frequent hoisting.

Intended Use

- Application of a hoist sling without rolling
- Minimise friction under pressure points
- Turning
- Positioning
- Move higher up in bed



Scan here to learn more!



LARGE 99437
SMALL 99IM437/5

ETAC 4WAY GLIDE LPL

A complete in bed system for turning and positioning of users with reduced mobility, high weight, pain or pressure issues.

Intended Use

- In bed system, user with reduced mobility
- In bed system, heavy user or user with pain or pressure ulcer
- Support horizontal transfer
- Can be used with support of a hoist
- Turning
- Positioning



Scan here to learn more!



99IM110200



ETAC MULTIGLIDE GLOVES

Versatile for a range of patient handling situations, Etac's MultiGlide Gloves are ideal for small re-positioning adjustments, movements or smoothing out bed sheets. Also useful when checking for adequate support when using pressure relieving cushions or mattresses, the glove facilitates placement of the hands between the user and mattress or cushion.



99-IM35



POSITIONING WEDGES

Etac Positioning Wedges are a multi-purpose product for use as a standalone positioning device or as an additional stability tool with other products. Ideal as a temporary side bolster to help prevent falls.

Key features

- Small and large sizes available
- Comfortably support patient in static position
- Minimises strain on care
- Non-slip cover for easy wipe cleaning



99-IM53NS



WONDER SHEET RANGE

The 'Wonder Sheet' fitted sheet is an effective therapeutic tool by increasing bed mobility. The design makes turning over and getting in and out of bed much easier. With your feet on the poly/cotton fabric at the bottom and your hips on the smooth satin fabric, it gives you the leverage to make turning over in bed and swinging legs in and out of bed much easier. This increased movement reduces pain and discomfort, and assists carers in moving the user.

Available in Single, Long Single, King Single, Double and Queen.
2 models available – Wonder Sheet and Wonder Sheet Plus.

Wondersheet Original



Wondersheet Plus



MANUAL HANDLING

SLIDE ASSIST



20-050

One Way Slide Tubes

Continuous loops that can be used on open or closed back chairs. Available in regular or bariatric sizes, for wheelchairs, for lounge or recliner tilt chairs, for tub or fall out chairs. For Velvet and Non Slip.

Small	37 x 43 cm Nonslip
Medium	44 x 45 cm Nonslip
Large	45 x 50 cm Nonslip
Bariatric	50 x 60 cm Nonslip 55 x 70 cm Nonslip 55 x 80 cm Nonslip



20-051

One Way Slide - Single Layer with Loop

Non-slip can be folded with their handle placed through the back of an open back chair, allowing the carer to use their strong leg muscles to easily position the person at the back of the chair or wheelchair.

Long single layer - sizing and model

Velour	130 x 40 cm
Nonslip	130 x 40 cm



20-0502

Easy Slide Gloves

One carer can move someone across a bed in three easy movements:

- 1/ Slide ankles across,
- 2/ Slide buttocks across,
- 3/ With patient's head on a pillow, slide shoulders across.

Two carers can easily slide someone up a bed, each with a glove on one hand. Place a glove under heels when sliding someone to reduce shear.

SLIDE SHEETS

Easy Slide sheets allows carers to slide the top sheet in any direction, or fold one sheet over to slide someone across a bed.



20-050_

Slide Sheet

Nylon, single-layer slide sheets are used to reposition a person in bed.

Available in 2 sizes

140 x 100 cm
140 x 200 cm



20-051_

One Way Slide Lite Velvet

Lite Velvet – a single layer One Way Slide used to help position a person on pressure cushions if they do not have severe decubitus problems or for positioning.




20-0508

Slide Sheet with Handles

Easy Slide Flat Sheet with Handles
Designed with additional Nylon webbing handgrips to assist the carer in gripping the sheet.

Learn more about other popular products that may support your home health care journey.



Statewide Home Health Care
909 Nepean Highway, Bentleigh 3204
T. (03) 9591 6234 | sales@shhc.com.au | shhc.com.au



Hire and sales also available on the Mornington Peninsula:

Peninsula Home Health Care
21 Sir Laurence Drive, Seaford 3198
T. 9786 7004 | sales@phhc.com.au | phhc.com.au

 @shhc_cares  @Statewide Home Health Care  @StatewideHomeHealthCare